

ФРЕЗЫ С КРУПНОЙ КРЕСТООБРАЗНОЙ НАРЕЗКОЙ



406102



406202



406302



406402



406502



406602



406001

Ø 1/10 mm	23	23	23	23	23	23	50
L mm	14,0	14,0	14,0	14,0	14,0	14,0	13,0
Тип нарезки	190	190	190	190	190	190	190
ISO 500 104	199 190 023	200 190 023	186 190 023	143 190 023	123 190 023	291 190 023	274 190 050

ФРЕЗЫ С КРУПНОЙ КРЕСТООБРАЗНОЙ НАРЕЗКОЙ



406801



406802



406901



407001



407101



407301



407401



407402

Ø 1/10 mm	40	40	50	60	60	60	60	50
L mm	14,0	12,0	14,0	13,0	12,0	14,5	14,0	13,0
Тип нарезки	190	190	190	190	190	190	190	190
ISO 500 104	201 190 040	201 190 040	200 190 050	274 190 060	237 190 060	277 190 060	201 190 060	143 190 050